

Annexure IPR Sheet	Doc. No.:	RSP/F/QA/24/05b-00
	Supersedes:	NIL
	Effective date:	24-12-2021
	Page 1 of 1	

PATIENT NAME: MR.KUSHAGRA TOKAS
DOCTOR: DR.GULVIR PONIYA


CASE ID: 1909071612
PLAN NO: 1

*	POSITIVE ATTACHMENT
—	NEGATIVE ATTACHMENT
▲	MISSING TOOTH
▨	IPR. (IN MM)
B	BITE RAMPS

@	BEFORE ALIGNER
●	PONTIC / BRIDGE / IMPLANT
X	EXTRACTION
◎	BUTTON ATTACHMENT

TOTAL NUMBER OF ALIGNER

STARTER KIT: YES



CATEGORY: AWD-25+

Remarks:

Attachments may vary during the course of treatment.

Note:-The Attachment Transfer template has a layer of isoflan sheet beneath it, Which is to be removed before using the Template.

Bite blocks (if present) are not to be transferred directly onto the tooth. They are incorporated within the aligners itself.

- Refinements may be required.
- Kindly ensure proper attachments transfer and their placements through-out the treatment.
- Use hand strips and/or Inter proximal reduction burs for Inter proximal reductions.
- Kindly ensure adequate patient compliance (minimum wear 20-22hrs./day)
- The Attachment Transfer Template has a layer of isoflan sheet beneath it, which has to be removed before using the template.

CONSULTANT ORTHODONTIST:
DATE:

Dr. Akshita Sharma
2019-09-28

APPROVED BY:
TIME:

Dr. Varun Grover
14:54:26